# INHEALTH

TECHNOLOGIES®

## Blom-Singer® Gel Cap Insertion System

### INSTRUCTIONS FOR USE

**Blom-Singer®**

voice restoration systems

37575-011 Effective January 2020 / 2020 年 1 月生效 / 2020 年 1 月発効 / Në fuqi nga janari 2020 / 2020 اعتباراً من شهر يناير / В сила от януари 2020 г. / Stupa na snagu u siječnju 2020. / V platnosti od ledna 2020 / Gældende fra januar 2020 / Met ingang van januari 2020 / Kehtib alates jaanuarist 2020 / Voimassa tammikuusta 2020 alkaen / Date d'entrée en vigueur : janvier 2020 / Gültig ab Januar 2020 / Σε ισχύ από τον Ιανουάριο του 2020 / בתוקף החל מינואר 2020 / Hatályos 2020. januártól / Gildir frá janúar 2020 / Valido da gennaio 2020 / 2020년 1월 발효 / Spēkā no 2020. gada janvāra / Įsigalioja 2020 m. sausio mėn. / Gjelder fra januar 2020 / Wchodzi w życie w styczniu 2020 r. / Vingente em janeiro de 2020 / Cu intrare în vigoare din ianuarie 2020 / Версия от января 2020 г. / Važi od januara 2020. / Platné od januára 2020 / Efectivo a partir de enero de 2020 / Gäller från januari 2020 / Geçerlilik Tarihi: Ocak 2020

Blom-Singer® and InHealth Technologies® are registered trademarks of Freudenberg Medical, LLC.

Patent(s): [www.inhealthpatents.com](http://www.inhealthpatents.com)



Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.



## INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

## DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagreame / Схеми / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar

## TABLE OF CONTENTS

目录 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /  
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /  
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /  
 Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins /  
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Gel Cap Insertion System	6	HU	Gélkapszula-behelyező rendszer	47
ZH	凝胶盖插入系统	8	IS	Ísetningarkerfi með hlauphyliki	50
JA	ゲルカプセル挿入システム	11	IT	Sistema di inserimento del cappuccio in gel	52
SQ	Sistemi i vendosjes të kapsulës xhelatinore	14	KO	젤라틴 캡슐 삽입 시스템	55
16	نظام إدخال كبسولة الجيلاتين	AR	LV	Gela vāciņu ievietošanas sistēma	58
BG	Система за поставяне на капачета от гел	19	LT	Gelinio dangtelio įterpimo sistema	60
HR	Sustav za umetanje s gel kapicom	22	NO	Innsettingssystem for gelkapsel	63
CS	Systém pro zavedení gelové krytky	24	PL	System wprowadzania kapturków żelowych	65
DA	Gelkapselindføringsssystem	27	PT	Sistema de Inserção de Cápsula de Gel	68
NL	Inbrengsysteem met gelkapje	29	RO	Sistem de introducere a capsulelor de gel	71
ET	Geelkapsli sisestussüsteem	32	RU	Система установки гелевого колпачка	73
FI	Geelisuojan asetusjärjestelmä	34	SR	Sistem za umetanje kapsula sa gelom	76
FR	Système d'insertion à capsule de gel	37	SK	Systém na nasunutie gélovej krytky	79
DE	Gelkappeneinführsystem	39	ES	Sistema de inserción por cápsula de gel	81
EL	Σύστημα εισαγωγής κάψουλας γέλης	42	SV	System för insättning av gelkapsel	84
45	ל'ג תולוספק תרדחהל תכרעמ	HE	TR	Jel Kapak Takma Sistemi	86

## ENGLISH

### **BLOM-SINGER® GEL CAP INSERTION SYSTEM** for Blom-Singer® Low Pressure Voice Protheses

#### **PRODUCT DESCRIPTION**

The Gel Cap Insertion System is designed to provide a method for inserting a 16 Fr. or 20 Fr. Blom-Singer Low Pressure Voice Prosthesis. Using the specially designed gel cap loading tool, the user places a capsule (gel cap) over the tip of the voice prosthesis to enable insertion into the tracheoesophageal puncture (opening made by surgeon between trachea and esophagus for placement of a voice prosthesis). The Gel Cap Insertion System (diagram 1) contains the following necessary items to insert a 16 Fr. or 20 Fr. Blom-Singer Low Pressure Voice Prosthesis: one (1) dual gel cap loading tool; ten (10) 16 Fr. gel caps; ten (10) 20 Fr. gel caps; one (1) 16 Fr. pushrod and one (1) 20 Fr. pushrod. Use only Blom-Singer 16 Fr. or 20 Fr. Gel Caps with this device. The gel cap loading tool and pushrods are reusable. The gel caps are for single use only; if the gel caps are used as intended, they dissolve and cannot be reused.

#### **INDICATIONS** (Reasons to prescribe the device or procedure)

The Blom-Singer Gel Cap Insertion System is designed for use by individuals who have had a laryngectomy and use a Blom-Singer 16 Fr. or 20 Fr. Low Pressure Voice Prosthesis.

#### **CONTRAINDICATIONS** (Reasons that make it inadvisable to prescribe the particular device or procedure)

The Blom-Singer Gel Cap Insertion System is a medical product that should be used only by individuals with experience and training in its use and care.

#### **WARNINGS AND PRECAUTIONS**

The Gel Cap Insertion System is a medically prescribed device. Instruction in its use must be provided by a qualified medical professional trained in the use of this particular device. The Blom-Singer Gel Cap Insertion System is designed for use with Blom-Singer Low Pressure Voice Protheses only. Do not use this product with any other type of device. The Gel Cap Insertion System does not function with special order low pressure voice protheses. Never use the pushrod as a voice prosthesis inserter. Only use the inserter provided with the prosthesis. Handle the voice prosthesis and gel cap loading tool carefully to avoid damage. Do not use a damaged voice prosthesis or gel cap loading tool. Use of a damaged device may result in unsuccessful insertion. If there are tears, cracks, or structural damage to the Gel Cap Insertion System or prosthesis, discontinue use. Failure to follow the instructions for use could result in difficulty or inability to insert or load the prosthesis.

#### **COMPLICATIONS**

Although rare, the following complications have been identified to occur with the use of the Blom-Singer Gel Cap Insertion System. They include:

- Tearing or other damage to the voice prosthesis from improper use of the gel cap loading tool.
- Aspiration of component(s).
- Allergic reaction to any of the device materials.

#### **INSTRUCTIONS FOR USE**

##### Loading the Gel Cap

*Please refer to the diagrams at the front of this instruction manual.*

The following procedure is provided by Eric D. Blom, Ph.D. as instructions for loading the voice prosthesis into the gel cap. The Gel Cap Insertion System, voice prosthesis, and user's hands should always be thoroughly cleaned to avoid introducing contaminants into the body. Gloves are

recommended. Be certain the tip of the voice prosthesis, gel cap loading tool, pushrod, and your hands are completely dry before beginning the following procedure. Failure to do so may cause the gel cap to dissolve prematurely.

Before beginning the following procedure, remove the items from the kit that correspond in size to your voice prosthesis, e.g., a 16 Fr. prosthesis uses the 16 Fr. gel caps and 16 Fr. pushrod.

1. Insert the neck strap of the voice prosthesis (diagram 2c) into the appropriate size hole of the dual gel cap loading tool on the side marked "top" (diagram 3). A 16 Fr. diameter prosthesis fits into the 16 Fr. **hexagon** hole and a 20 Fr. diameter prosthesis fits into the 20 Fr. **round** hole.
2. Holding the neck strap, gently pull the prosthesis down and through the hole (diagram 4) until the esophageal flange is positioned over the hole (diagram 5).
3. Grasping the body of the voice prosthesis at the base of the neck strap (diagram 2b), slowly and gently pull the prosthesis (diagram 6) until the esophageal flange on the voice prosthesis begins to fold forward inside the hole (diagram 7). Do not pull the prosthesis beyond this point. If the prosthesis is over-pulled and releases from the gel cap loading tool, return to step 1.
4. From the gel cap container, remove the 16 Fr. or 20 Fr. gel cap that corresponds to the diameter of your voice prosthesis. Discard the longer end of the capsule (diagram 8). Place the open end of the gel cap in the hole of the loading tool with the folded flange of the prosthesis (diagram 9). The gel cap should fit into the groove and cover the prosthesis tip. Verify that the gel cap is securely in place.
5. Fully insert the appropriate-sized pushrod into the center of voice prosthesis at the base of the neck strap. Place a finger on the tip of the gel cap to hold it in place while simultaneously pushing the voice prosthesis from the bottom (diagram 10). Push gently until the folded esophageal flange is fully positioned in the gel cap.
6. Remove your fingertip from the gel cap and push gently with the pushrod in the direction indicated by the arrow in diagram 11.
7. Continue to push the prosthesis until it is completely out of the gel cap loading tool (diagram 12). Hold the prosthesis at the base of the neck strap and remove the pushrod from the prosthesis. The voice prosthesis is now fully loaded in the gel cap and ready for insertion in the usual manner. Refer to the Instruction for Use Manual provided with your voice prosthesis for device insertion instructions.

#### **CLEANING AND CARE INSTRUCTIONS**

Wash the dual gel cap loading tool and the pushrod thoroughly with a mild detergent and rinse with water. Dry the devices with gauze or lint-free cloth. Never use facial or toilet tissue to dry the dual gel cap loading tool, inserter, or voice prosthesis, as particles of lint or fabric could be inhaled into the airway upon insertion of the prosthesis. Once the dual gel cap loading tool and the pushrod have been cleaned and thoroughly dried, they may be stored in a clean, resealable plastic bag.

#### **SPECIAL STORAGE AND/OR HANDLING CONDITIONS**

**Only the gel caps packaged with the device have special storage and/or handling conditions.** They must be **kept dry** and **out of direct sunlight**.

#### **DISPOSAL INSTRUCTIONS**

This product is not biodegradable and must be considered contaminated upon use. Carefully dispose of this device per local guidelines.

## ORDERING INFORMATION

USA

Online Ordering [www.inhealth.com](http://www.inhealth.com)

USA

InHealth Technologies

1110 Mark Avenue

Carpinteria, CA 93013-2918

Customer Service: +1-800-477-5969

Fax: +1-888-371-1530

Email: [order@inhealth.com](mailto:order@inhealth.com)

INTERNATIONAL

Contact InHealth Technologies Customer Service at +1-805-684-9337 for international distributor contact.

## PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact

**[productcomplaints@inhealth.com](mailto:productcomplaints@inhealth.com)**.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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### SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Batch code; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified

**SYMBOLS GLOSSARY (USA only)**

<b>SYMBOL</b>	<b>STANDARD No.</b>	<b>SYMBOL TITLE &amp; REF No.</b>	<b>MEANING OF SYMBOL</b>
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaged has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
	EN ISO 15223-1:2016	Keep away from sunlight; 5.3.2	Indicates a medical device that needs protection from light sources
	EN ISO 15223-1:2016	Keep dry; 5.3.4	Indicates a medical device that needs to be protected from moisture
	ISO 7000	Packaging unit; 2794	To indicate the number of pieces in the package

# INHEALTH

## TECHNOLOGIES®

---

### UNITED STATES

1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
Toll-Free: 800.477.5969  
Toll-Free Fax: 888.371.1530

### INTERNATIONAL

Tel: +1.805.684.9337  
Fax: +1.805.684.8594  
E-mail: [order@inhealth.com](mailto:order@inhealth.com)

### ORDER ONLINE

[www.inhealth.com](http://www.inhealth.com)

**Freudenberg Medical, LLC**  
1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
[www.inhealth.com](http://www.inhealth.com)

**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands